

## COVID-19 Ag TEST



### INTENDED USE

Turklab's COVID-19 Ag Test is a rapid immunochromatographic assay for the qualitative detection of SARS-CoV-2 nucleocapsid antigen in human nasal swabs in a professional environment.

- Reliable results in the acute phase (first 5-7 days) of COVID-19 infection.
- Varified detection of Omicron variant's nucleocapsid antigen.
- Showed no reactivity with other human coronavirus species or other respiratory viruses.

### PERFORMANCE EVALUATION

- Sensitivity : 96,58%
- Specificity: 99,47%



PRODUCT	TEST TYPE	SAMPLE TYPE	READING TIME	KIT SIZE	KIT COMPONENTS
COVID-19 Ag	CASSETTE	NASAL SWAB	20 Min.	25	Test Cassette + Sterile Swab + Sample Extraction Reagent Tube + Instructions For Use